

# CTA-100 | DIGITAL EYEPIECE ADAPTER



**Spend less money to get a Digital Slit Lamp.**  
**Digital eyepiece adapters take optical equipment into a new digital era.**



DIA: 38MM

Small Eyepiece Connector is Available now, Dedicated to Slit Lamp with smaller diameter eyepieces, such as Topcon, Haag Streit.



Adjustable support case & Camera position could hold the size of Smart-phones from 4" to 5.5" .



● Slit Lamp



● Telescope



● Microscope

CTA-100 Eyepiece Adapter converts smart phones into professional digital eyepieces for the traditional optical Instruments. Just attach the equipped smart phone to the slit lamps, microscopes, telescopes, Spotting Scope, Binocular or Monocular etc. With built-in camera app, you can take sharp images, capture videos, and then share all information via email, SMS, twitter, facebook.

### Features:

1. Adjustable support case dimension makes it available for the most mainstream smart-phones. The sizes of the smart-phones can be from 4" to 5.5" .
2. Ultra-durable ABS plastic.
3. Fast Installation, Easy Operation
4. Send out high-res images directly from your phone via email, SMS, twitter or 3rd party apps.
5. It is designed to work with eyepieces that are 1in – 1.8in in diameter (25mm – 48mm).
6. Silicone sizing bands are included in the case to fit for the smaller diameter eye-pieces.
7. Work with

- Slit Lamp
- Microscopes
- Telescopes
- Binoculars
- Spotting scopes
- Other Optical Instruments.

### ● Installation Diagram



### ● Standard Accessories:



Big Eyepiece Connector: φ 48mm      Small Eyepiece Connector: φ 36mm



Big Silicone Sizing bands:1PC      Small Silicone Sizing bands:2PCS