



Daud Jee Mfg. Co.

PLASTIC SURGERY

INSTRUMENTS





NEEDLE HOLDERS

Micron Fine Needle Holder

- . 12.5cm/5" length
- . Diamond dusted jaws



DJE-1010



Derf Needle Holder

- . 12cm/4¾" length
- . Tungsten carbide jaws



DJE-1011

Smooth



DJE-1012



Halsey Needle Holder

- . 13cm/5" length
- . Tungsten carbide jaws



DJE-1013

Smooth



DJE-1014



Crile Wood Needle Holder

- . Tungsten carbide jaws



DJE-1015
15cm/6"



DJE-1016
18cm/7"



DJE-1017
20cm/8"






NEEDLE HOLDERS

Nieverts Needle Holder (offset ring)

- . 13cm/5¼" length
- . Tungsten carbide jaws



DJE-1018

Foster Gillies Combined Needle Holder/Scissors

- . 13cm/5" length
- . Tungsten carbide jaws



DJE-1019

Gillies Combined Needle Holder/Scissors

- . 16cm/6¼" length
- . Tungsten carbide jaws



DJE-1020



NEEDLE HOLDERS

Webster

DJE-1021 11.5 cm, smooth
DJE-1022 13.0 cm, smooth, TC



Halsey

DJE-1023 13.0 cm, smooth, TC
DJE-1024 13.0 cm, cross serrated, TC



Par

DJE-1025 11.5 cm, cross serrated, TC
DJE-1026 13.0 cm, cross serrated, TC





NEEDLE HOLDERS

Derf

DJE-1027 12.0 cm, serrated
DJE-1028 12.0 cm, cross serrated, TC



Neivert

DJE-1029 13.0 cm, cross serrated, TC



Swedish

DJE-1030 15.5 cm, cross serrated, TC






NEEDLE HOLDERS

Ryder



- DJE-1031 13.0 cm, cross serrated, TC
- DJE-1032 15.0 cm, cross serrated, TC
- DJE-1033 18.0 cm, cross serrated, TC

Micro Ryder



- DJE-1034 13.0 cm, cross serrated, TC
- DJE-1035 15.0 cm, cross serrated, TC
- DJE-1036 18.0 cm, cross serrated, TC

Euphrate-Pasque

. extremely delicate



- DJE-1037 14.0 cm, cross serrated, TC
- DJE-1038 15.0 cm, cross serrated, TC
- DJE-1039 18.0 cm, cross serrated, TC
- DJE-1040 20.0 cm, cross serrated, TC




NEEDLE HOLDERS

Micro-Vascular




- DJE-1041 15.0 cm, cross serrated, TC
- DJE-1042 18.0 cm, cross serrated, TC

Crile-Wood



- DJE-1043 15.0 cm, smooth, TC
- DJE-1044 15.0 cm, cross serrated, TC
- DJE-1045 18.0 cm, cross serrated, TC
- DJE-1046 20.0 cm, cross serrated, TC
- DJE-1047 23.0 cm, cross serrated, TC

Lemmon




- DJE-1048 14.0 cm, cross serrated, TC
- DJE-1049 17.0 cm, cross serrated, TC
- DJE-1050 19.0 cm, cross serrated, TC



NEEDLE HOLDERS

Baumgartner




DJE-1051 12.5 cm, cross serrated, TC
DJE-1052 14.0 cm, cross serrated, TC

Mayo Hegar



DJE-1053 14.0 cm, cross serrated, TC
DJE-1054 16.0 cm, cross serrated, TC
DJE-1055 18.0 cm, cross serrated, TC
DJE-1056 20.0 cm, cross serrated, TC

Wire Twister



DJE-1057 13.0 cm, round jaw, cross serrated, TC
DJE-1058 16.0 cm, blunt jaw, cross serrated, TC
DJE-1059 18.0 cm, round jaw, cross serrated, TC
DJE-1060 18.0 cm, blunt jaw, cross serrated, TC



NEEDLE HOLDERS

Corwin Wire Twister



- DJE-1061 16.0 cm, long round jaw, cross serrated, TC
- DJE-1062 16.0 cm, blunt jaw, cross serrated, TC
- DJE-1063 16.0 cm, round jaw, cross serrated, TC

Olsen Hegar

. needleholder/scissors combination



- DJE-1064 14.0 cm, cross serrated, TC
- DJE-1065 18.0 cm, cross serrated, TC

Gillies

. needleholder/scissors combination



- DJE-1066 16.0 cm, cross serrated, curved, TC



NEEDLE HOLDERS

Barraquer Micro

. with round handle, without lock, tip 0.8 × 1.2 mm

DJE-1067	14.0 cm, straight, diamond dusted
DJE-1068	14.0 cm, curved, diamond dusted
DJE-1069	15.0 cm, straight, diamond dusted
DJE-1070	15.0 cm, curved, diamond dusted
DJE-1071	18.0 cm, straight, diamond dusted
DJE-1072	18.0 cm, curved, diamond dusted



Barraquer Micro

. with round handle, with lock, tip 0.8 × 1.2 mm

DJE-1073	14.0 cm, straight, diamond dusted
DJE-1074	14.0 cm, curved, diamond dusted
DJE-1075	15.0 cm, straight, diamond dusted
DJE-1076	15.0 cm, curved, diamond dusted
DJE-1077	18.0 cm, straight, diamond dusted
DJE-1078	18.0 cm, curved, diamond dusted





NEEDLE HOLDERS

Barraquer Delicate

. with round handle, without lock, tip 1.2 × 1.8 mm

- DJE-1079 14.0 cm,
cross serrated, straight, TC
- DJE-1080 14.0 cm,
cross serrated, curved, TC
- DJE-1081 15.0 cm,
cross serrated, straight, TC
- DJE-1082 15.0 cm,
cross serrated, curved, TC
- DJE-1083 18.0 cm,
cross serrated, straight, TC
- DJE-1084 18.0 cm,
cross serrated, curved, TC



Barraquer Micro

. with round handle, without lock, tip 0.8 × 1.2 mm

- DJE-1085 14.0 cm,
cross serrated, straight, TC
- DJE-1086 14.0 cm,
cross serrated, curved, TC
- DJE-1087 15.0 cm,
cross serrated, straight, TC
- DJE-1089 15.0 cm,
cross serrated, curved, TC
- DJE-1090 18.0 cm,
cross serrated, straight, TC
- DJE-1091 18.0 cm,
cross serrated, curved, TC





NEEDLE HOLDERS

Barraquer Delicate

. with round handle, with lock, tip 1.2 × 1.8 mm



DJE-1092	14.0 cm, cross serrated, straight, TC
DJE-1093	14.0 cm, cross serrated, curved, TC
DJE-1094	15.0 cm, cross serrated, straight, TC
DJE-1095	15.0 cm, cross serrated, curved, TC
DJE-1096	18.0 cm, cross serrated, straight, TC
DJE-1097	18.0 cm, cross serrated, curved, TC

Barraquer Micro

. with round handle, with lock, tip 0.8 × 1.2 mm



DJE-1098	14.0 cm, cross serrated, straight, TC
DJE-1099	14.0 cm, cross serrated, curved, TC
DJE-1100	15.0 cm, cross serrated, straight, TC
DJE-1101	15.0 cm, cross serrated, curved, TC
DJE-1102	18.0 cm, cross serrated, straight, TC
DJE-1103	18.0 cm, cross serrated, curved, TC



NEEDLE HOLDERS

Castroviejo

. with flat handle, with lock



DJE-1104	14.0 cm, smooth, straight, TC
DJE-1105	14.0 cm, cross serrated, straight, TC
DJE-1106	14.0 cm, smooth, curved, TC
DJE-1107	14.0 cm, cross serrated, curved, TC
DJE-1108	18.0 cm, smooth, straight, TC
DJE-1109	18.0 cm, cross serrated, straight, TC
DJE-1110	18.0 cm, smooth, curved, TC
DJE-1111	18.0 cm, cross serrated, curved, TC

Castroviejo

. with flat handle, without lock



DJE-1112	14.0 cm, straight, diamond dusted
DJE-1113	14.0 cm, curved, diamond dusted
DJE-1113	18.0 cm, straight, diamond dusted
DJE-1114	18.0 cm, curved, diamond dusted



NEEDLE HOLDERS

Castroviejo

. with flat handle, with lock



- DJE-1115 14.0 cm, straight, diamond dusted
- DJE-1116 14.0 cm, curved, diamond dusted
- DJE-1117 18.0 cm, straight, diamond dusted
- DJE-1118 18.0 cm, curved, diamond dusted

Castroviejo Ryder

. with flat handle, with lock



- DJE-1119 14.0 cm, cross serrated, straight, TC
- DJE-1120 18.0 cm, cross serrated, straight, TC

Barraquer

. with round handle, with inserted jaws 0.75mm



- DJE-1121 12.0 cm, straight, without lock, smooth
- DJE-1122 12.0 cm, curved, without lock, smooth



FACE LIFT SCISSORS

Gorney Face Lift Scissors

- . Tungsten carbide edges

Straight	Curved	Length
DJE-1123	DJE-1124	17cm/6¾"
DJE-1125	DJE-1126	19cm/7½"
DJE-1127	DJE-1128	23cm/9"



Flat tips



Gorney-Freeman Face Lift Scissors

- . Tungsten carbide edges
- . Open shanks

Straight	Curved	Length
DJE-1129	DJE-1130	17cm/6¾"
DJE-1131	DJE-1132	19cm/7½"
DJE-1133	DJE-1134	23cm/9"



Flat tips



Kaye Face Lift Scissors

- . Tungsten carbide edges
- . Curved

Straight	Curved
DJE-1135	17cm/6¾"
DJE-1136	19cm/7½"
DJE-1137	23cm/9"



Aufricht Scissors

- . 14cm/5¼" length
- . Tungsten carbide edges
- . Curved

DJE-1138



Goldman-Fox Scissors

- . 13cm/5" length
- . Tungsten carbide edges
- . Tapered blades to fine tips
- . Curved

DJE-1139





DISSECTING SCISSORS

Ragnell (Kilner) Scissors

- . Tungsten carbide edges
- . Curved flat tips

Code	Length
DJE-1140	13cm/5¼"
DJE-1141	15cm/6"
DJE-1142	18cm/7"



Greenberg PAR Scissors

- . 11,5cm/4¼" length
- . Tungsten carbide edges
- . Curved tapered flat tips

DJE-1143



Littler Plastic Surgery Scissors

- . 12cm/4¾" length
- . Curved
- . Hole in the blades for the suture to be threaded through

DJE-1144



Kahn Face Lift Scissors

- . 14,5cm/5¾" length
- . Super Cut edges
- . Curved

DJE-1145



Fomon Saber Face Lift Scissors

- . 13,5cm/5¼" length
- . Super Cut edges
- . Curved, round blunt tips

DJE-1146





SCISSORS

Iris Scissors

- . 11.5cm/4½" length
- . Straight
- . Micro serrated cutting edge

DJE-1147



Iris Scissors

- . 11.5cm/4½" length
- . Curved
- . Micro serrated cutting edge

DJE-1148



Strabismus Scissors

- . 11.5cm/4½" length
- . Straight
- . Micro serrated cutting edge

DJE-1149



Strabismus Scissors

- . 11.5cm/4½" length
- . Curved
- . Micro serrated cutting edge

DJE-1150



Waldram Scissors (with larger bows)

- . 12.5cm/5" length
- . Curved
- . Micro serrated cutting edge

DJE-1151





SCISSORS

Gorney



- DJE-1152 18.0 cm, curved, serrated, TC
- DJE-1153 20.0 cm, straight, serrated, TC
- DJE-1154 20.0 cm, curved, serrated, TC
- DJE-1155 23.0 cm, curved, serrated, TC

Gorney



- DJE-1156 20.0 cm, straight, serrated
- DJE-1157 20.0 cm, curved, serrated

Freeman Gorney



- DJE-1158 20.0 cm, straight, TC
- DJE-1159 20.0 cm, curved, TC

Freeman



- DJE-1160 18.0 cm, straight, TC
- DJE-1161 18.0 cm, curved, TC



SCISSORS

Aston-Face



DJE-1162 17.0 cm, curved
DJE-1163 20.0 cm, curved

Rees



DJE-1163 17.0 cm, straight, serrated, TC
DJE-1164 17.0 cm, curved, serrated, TC
DJE-1165 20.0 cm, straight, serrated, TC
DJE-1166 20.0 cm, curved, serrated, TC

Kaye



DJE-1167 15.0 cm, curved, serrated, TC
DJE-1168 23.0 cm, curved, serrated, TC

Kaye



DJE-1169 11.5 cm, curved, serrated, TC



SCISSORS

Freeman (Kaye)



DJE-1170 18.0 cm, curved, serrated, TC
DJE-1171 20.0 cm, curved, serrated, TC

Davis



DJE-1172 18.0 cm, curved,
serrated, sharp/sharp, TC

Castanares



DJE-1173 16.5 cm, straight, TC

Goldman Fox



DJE-1174 13.0 cm, curved, serrated, TC



SCISSORS

Par Style

. very delicate



DJE-1175 11.5 cm, curved, TC
DJE-1176 14.5 cm, curved, TC
for serrated blades please add "WS"

Stevens Tenotomy



DJE-1177 11.5 cm, straight, TC
DJE-1178 11.5 cm, curved, TC

Stevens Tenotomy



DJE-1179 11.0 cm, straight, sharp/sharp
DJE-1180 11.0 cm, curved, sharp/sharp
DJE-1181 11.0 cm, straight, blunt/blunt
DJE-1182 11.0 cm, curved, blunt/blunt

Devine



DJE-1183 fine dissecting scissors
with sharp tips,
11.0 cm, curved, TC



SCISSORS

Kilner

- DJE-1184 12.0 cm, straight, TC
- DJE-1185 12.0 cm, curved, TC
- DJE-1186 12.0 cm, curved
1 blade TC
1 blade serrated



Strabismus

- DJE-1187 11.0 cm, straight, TC
- DJE-1188 11.0 cm, curved, TC
- DJE-1189 11.0 cm, straight, supercut
- DJE-1190 11.0 cm, curved, supercut
for serrated blades please add "WS"



Iris

- DJE-1191 11.5 cm, round, straight, sharp/sharp
- DJE-1192 11.5 cm, round, curved, sharp/sharp
- DJE-1193 11.5 cm, round, straight, sharp/sharp, TC
- DJE-1194 11.5 cm, round, curved, sharp/sharp, TC
for serrated blades please add "WS"



Iris

- DJE-1195 11.5 cm, square, straight, sharp/sharp
- DJE-1196 11.5 cm, square, curved, sharp/sharp
- DJE-1197 11.5 cm, square, straight, sharp/sharp, TC
- DJE-1198 11.5 cm, square, curved, sharp/sharp, TC
for serrated blades please add "WS"





SCISSORS

Joseph



- DJE-1199 14.0 cm, straight, sharp/sharp
- DJE-1200 14.0 cm, curved, sharp/sharp
- DJE-1201 14.0 cm, straight, sharp/sharp, TC
- DJE-1202 14.0 cm, curved, sharp/sharp, TC

Locklin



- DJE-1203 16.0 cm, straight, angled, sharp/sharp, TC

Shea



- DJE-1204 12.0 cm, sharp, fine, curved



SCISSORS

Heymann



DJE-1205 18.0 cm, angled to side

Knight



DJE-1206 18.0 cm, angled to side, serrated,
delicate blades, TC

Wilmer Converse



DJE-1207 10.5 cm, angled up



SCISSORS

Fomon



DJE-1208 13.5 cm, lower lateral scissor, strongly curved

Fomon



DJE-1209 14.5 cm, angled, TC

Cottle



DJE-1210 17.0 cm, angled, TC



SCISSORS

Caplan Septum



- DJE-1211 22.0 cm, one blade serrated, double action, straight
- DJE-1212 22.0 cm, one blade serrated, double action, curved

Metzenbaum blunt/blunt

- DJE-1213 14.5 cm, straight, TC
 - DJE-1214 14.5 cm, curved, TC
 - DJE-1215 18.0 cm, straight, TC
 - DJE-1216 18.0 cm, curved, TC
 - DJE-1217 20.5 cm, straight, TC
 - DJE-1218 20.5 cm, curved, TC
 - DJE-1219 23.0 cm, straight, TC
 - DJE-1220 23.0 cm, curved, TC
 - DJE-1221 28.0 cm, straight, TC
 - DJE-1222 28.0 cm, curved, TC
- for serrated blades please add "WS"



Metzenbaum, blunt/blunt

- DJE-1223 14.5 cm, straight, supercut
- DJE-1224 14.5 cm, curved, supercut
- DJE-1225 18.0 cm, straight, supercut
- DJE-1226 18.0 cm, curved, supercut
- DJE-1227 23.0 cm, straight, supercut
- DJE-1228 23.0 cm, curved, supercut





SCISSORS

Metzenbaum Fino blunt/blunt

DJE-1229	14.5 cm, straight, TC
DJE-1230	14.5 cm, curved, TC
DJE-1231	18.0 cm, straight, TC
DJE-1232	18.0 cm, curved, TC
DJE-1233	20.0 cm, straight, TC
DJE-1234	20.0 cm, curved, TC
DJE-1235	23.0 cm, straight, TC
DJE-1236	23.0 cm, curved, TC
DJE-1237	28.0 cm, straight, TC
DJE-1238	28.0 cm, curved, TC

for serrated blades please add "WS"



Metzenbaum sharp/sharp

DJE-1239	14.5 cm, straight, TC
DJE-1240	14.5 cm, curved, TC
DJE-1241	18.0 cm, straight, TC
DJE-1242	18.0 cm, curved, TC
DJE-1243	23.0 cm, straight, TC
DJE-1244	23.0 cm, curved, TC

for serrated blades please add "WS"



Mayo-Stille

DJE-1245	15.0 cm, straight, TC
DJE-1246	15.0 cm, curved, TC
DJE-1247	17.0 cm, straight, TC
DJE-1248	17.0 cm, curved, TC

for serrated blades please add "WS"





SCISSORS

Mayo-Stille



DJE-1249 17.0 cm, straight, supercut
DJE-1250 17.0 cm, curved, supercut

Mayo-Lexer



DJE-1251 21.0 cm, straight, TC
DJE-1252 21.0 cm, curved, TC

Operating sharp/blunt



DJE-1253 14.5 cm, straight, round, TC
DJE-1254 14.5 cm, curved, round, TC



SCISSORS

Operating sharp/sharp



DJE-1255 14.5 cm, straight, round, TC
DJE-1256 14.5 cm, curved, round, TC

Operating blunt/blunt



DJE-1257 14.5 cm, straight, round, TC
DJE-1258 14.5 cm, curved, round, TC

Universal Wire Scissors



DJE-1259 12.0 cm, serrated, TC

Stevens Tenotomy Scissors



DJE-1260 9.5 cm, curved, sharp/sharp, ribbon type
DJE-1261 12.0 cm, curved, sharp/sharp, ribbon type



SCISSORS

Iris

- DJE-1262 9.5 cm, straight, sharp/sharp, ribbon type
- DJE-1263 09.5 cm, curved, sharp/sharp, ribbon type
- DJE-1264 10.5 cm, straight, sharp/sharp, ribbon type
- DJE-1265 10.5 cm, curved, sharp/sharp, ribbon type



Littler

- DJE-1266 12.0 cm, curved, tips with eye for suture
- DJE-1267 14.0 cm, curved, tips with eye for suture



Spencer

- DJE-1268 9.0 cm, straight, ligature scissor



Knowles Bandage

- DJE-1269 14.0 cm, angled, serrated





SCISSORS

Lister Bandage



DJE-1270 14.0 cm, angled
DJE-1271 18.0 cm, angled

Freeman Kaye



DJE-1272 17.0 cm, curved

Micro Scissors

. spring scissors with round handle, blade 14 mm



DJE-1273 15.0 cm, straight
DJE-1274 15.0 cm, straight, blades serrated
DJE-1275 15.0 cm, curved
DJE-1276 15.0 cm, curved, blades serrated
DJE-1277 18.0 cm, straight
DJE-1278 18.0 cm, straight, blades serrated
DJE-1279 18.0 cm, curved
DJE-1280 18.0 cm, curved, blades serrated



SCISSORS

Slim Baby Metzenbaum Easy-Cut Scissors

- . 14cm/5½" length
- . Curved
- . Micro serrated cutting edge



DJE-1281



Slim Metzenbaum Easy-Cut Scissors

- . 18cm/7" length
- . Curved
- . Micro serrated cutting edge



DJE-1282



Baby Metzenbaum Easy-Cut Scissors

- . 14cm/5½" length
- . Curved
- . Micro serrated cutting edge



DJE-1283



McIndoe / Metzenbaum Easy-Cut Scissors

- . 18cm/7" length
- . Curved
- . Micro serrated cutting edge



DJE-1284



Crayford-Reynolds Scissors

- . 18cm/7" length
- . Curved
- . Micro serrated cutting edge



DJE-1285




SCISSORS

Mayo Easy-Cut Scissors

- . 14cm/5½" length
- . Straight
- . Micro serrated cutting edge



DJE-1286




Mayo Easy-Cut Scissors

- . 14cm/5½" length
- . Curved
- . Micro serrated cutting edge



DJE-1287




Mayo Easy-Cut Scissors

- . 17cm/6¾" length
- . Straight
- . Micro serrated cutting edge



DJE-1288



Mayo Easy-Cut Scissors

- . 17cm/6¾" length
- . Curved
- . Micro serrated cutting edge



DJE-1289





TENOTOMY SCISSORS

Stevens Tenotomy Scissors

- . 11.5cm/4½" length
- . Curved
- . Micro serrated cutting edge

DJE-1290



Stevens-Stille Pattern Tenotomy Scissors

- . 12.5cm/5" length
- . Curved
- . Micro serrated cutting edge
- . Larger oversized rings

DJE-1291



Jameson-Werber Scissors

- . 13.5cm/5¼" length
- . Curved
- . Micro serrated cutting edge

DJE-1292



Jameson Scissors

- . 15cm/6" length
- . Curved
- . Micro serrated cutting edge

DJE-1293



Jameson Scissors

- . 15cm/6" length
- . Curved
- . Strong pattern
- . Micro serrated cutting edge

DJE-1294






DISSECTING FORCEPS

Gillies Dissecting Forceps


- . 15cm/6" length
- . Toothed



DJE-1295

McIndoes Dissecting Forceps

- . 15cm/6" length
- . Non toothed



DJE-1296

Debakey-Diethrich Tissue Forceps


- . 15cm/6" length
- . 1.5mm tips



DJE-1297

Debakey Tissue Forceps

- . 15cm/6" length
- . 2mm tips



DJE-1298



DISSECTING FORCEPS

Bipolar Forceps

- . 12.5cm/5" length
- . 1mm tips



DJE-1299

Bipolar Forceps

- . 15cm/6" length
- . 1.5mm tips



DJE-1300

Bipolar Cable

- . 3m cable length
- . Block end/Banana plug fitting to machine end



DJE-1301

Adson Forceps

- . 12cm/4¾" length
- . Non toothed



DJE-1302



Adson Forceps

- . 12cm/4¾" length
- . Toothed



DJE-1303





DISSECTING FORCEPS

Gillies Dissecting Forceps

- . 15cm/6" length
- . Toothed



DJE-1304



McIndoes Dissecting Forceps

- . 15cm/6" length
- . Non toothed



DJE-1305



Debakey-Diethrich Tissue Forceps

- . 15cm/6" length
- . 1.5mm tips




DJE-1306



Debakey Tissue Forceps

- . 15cm/6" length
- . 2mm tips



DJE-1307





FACELIFT FORCEPS

Assumpcao Face Lift Forceps

- . 16cm/6¼" length
- . For precise marking of dissected facial flaps for excision



Without ratchet
DJE-1308

With ratchet
* DJE-1309

illustrated with ratchet *

Pitanguy Forceps

- . Flap demarcator for accurate marking of facial flaps for excision



Code	Length		
* DJE-1310	15cm/6"		
DJE-1311	26cm/10¼"	(for abdominal flaps)	illustrated*

Marten Flap Forceps

- . 14.5cm/5¾" length
- . For accurate marking of skin flaps



With ratchet
DJE-1312

Without ratchet
* DJE-1313

illustrated *

Martel Flap Grasping Forceps

- . 14cm/5½" length
- . Two pins in opposing holes on the jaws



DJE-1314



FORCEPS

Hartmann Mosquito

DJE-1315 10.0 cm, straight
DJE-1316 10.0 cm, curved



Halstead Mosquito

DJE-1317 12.5 cm, straight, 1x2 teeth
DJE-1318 12.5 cm, curved, 1x2 teeth
DJE-1319 12.5 cm, straight
DJE-1320 12.5 cm, curved



Kelly

DJE-1321 14.0 cm, straight
DJE-1322 14.0 cm, curved



Crile

. hemostatic forceps

DJE-1323 14.0 cm, straight
DJE-1324 14.0 cm, curved
DJE-1325 16.0 cm, curved





FORCEPS

Rochester Pean

. hemostatic forceps

DJE-1326	14.0 cm, straight
DJE-1327	16.0 cm, straight
DJE-1328	20.0 cm, straight
DJE-1329	16.0 cm, curved
DJE-1330	18.0 cm, curved
DJE-1331	20.0 cm, curved



Rochester-Ochsner

. hemostatic forceps

DJE-1332	14.0 cm, straight, 1x2 teeth
DJE-1333	14.0 cm, curved, 1x2 teeth
DJE-1334	16.0 cm, straight, 1x2 teeth
DJE-1335	16.0 cm, curved, 1x2 teeth



Adson

. hemostatic forceps

DJE-1336	18.5 cm, straight
DJE-1337	18.5 cm, curved





FORCEPS

Lahey Vulsellum



DJE-1338 15.0 cm
DJE-1339 20.0 cm

Allis



DJE-1340 14.0 cm, 4:5 teeth, Allis Baby
DJE-1341 15.0 cm, 4:5 teeth
DJE-1342 15.0 cm, 5:6 teeth
DJE-1343 19.0 cm, 5:6 teeth

Backhaus Towel



DJE-1344 9.0 cm
DJE-1345 11.0 cm
DJE-1346 13.0 cm
DJE-1347 15.0 cm



FORCEPS

Lorna Towel



DJE-1348 11.0 cm
DJE-1349 13.0 cm

Jones Towel



DJE-1350 9.0 cm

Foerster Sponge Holding



DJE-1351 18.0 cm
DJE-1352 20.0 cm
DJE-1353 25.0 cm



FORCEPS / THUMB

Adson



DJE-1354 12.0 cm, cross serrated, TC
DJE-1355 15.0 cm, cross serrated, TC

Adson



DJE-1356 12.0 cm, serrated, Titanium

Adson



DJE-1357 12.0 cm, serrated
DJE-1358 12.0 cm, serrated, Micro

Adson Dressing



DJE-1359 12.0 cm, cross serrated



FORCEPS / THUMB

Adson



DJE-1360 12.0 cm, cross serrated, 1x2 teeth, TC
DJE-1361 15.0 cm, cross serrated, 1x2 teeth, TC

Adson



DJE-1362 12.0 cm, smooth, 1x2 teeth, Titanium

Adson



DJE-1363 12.0 cm, smooth, 1x2 teeth

Adson Tubiana



DJE-1364 12.0 cm, 1x2 teeth, with platform 1 mm



FORCEPS / THUMB

Brown-Adson



DJE-1365 12.0 cm, cross serrated, Brown teeth, TC
DJE-1366 15.0 cm, cross serrated, Brown teeth, TC

Brown-Adson



DJE-1367 12.0 cm, Brown teeth

Adson Micro



DJE-1368 12.0 cm, cross serrated, 1x2 teeth, TC
DJE-1369 15.0 cm, cross serrated, 1x2 teeth, TC

Par



DJE-1370 12.0 cm, cross serrated, TC



FORCEPS / THUMB

Adson Delicate



DJE-1371 12.0 cm, smooth, TC
DJE-1372 15.0 cm, cross serrated, TC

Gorney-Adson



DJE-1373 12.0 cm, 1x2 teeth,
with platform, diamond dusted

Castroviejo



DJE-1374 10,5 cm, 1.5 mm tip, 1x2 teeth, with platform
DJE-1375 10,5 cm, 3 mm tip, 1x2 teeth, with platform
DJE-1376 10,5 cm, 5 mm tip, 1x2 teeth, with platform
DJE-1377 10,5 cm, 9 mm tip, 1x2 teeth, with platform

Castroviejo-Brown



DJE-1378 11.0 cm, brown teeth



FORCEPS / THUMB

Graefe Iris



DJE-1379 10.0 cm, smooth, 1x2 teeth

Graefe Iris



DJE-1380 10.0 cm, serrated

Julian Suture



DJE-1381 11.0 cm, smooth, curved

Adson-Graefe



DJE-1382 12.0 cm, multi-toothed



FORCEPS / THUMB

Bishop-Harmon



- DJE-1383 8.5 cm, smooth, 1x2 teeth, fine
- DJE-1384 8.5 cm, smooth, 1x2 teeth

Green Fixation

- DJE-1385 10.0 cm, 5 mm tip, without lock
- DJE-1386 10.0 cm, 10 mm tip, without lock
- DJE-1387 10.0 cm, 12 mm tip, without lock
- DJE-1388 10.0 cm, 5 mm tip, with lock
- DJE-1389 10.0 cm, 10 mm tip, with lock
- DJE-1390 10.0 cm, 12 mm tip, with lock



Semken



- DJE-1391 15.0 cm, cross-serrated, 1x2 teeth, TC

McIndoe Dissecting



- DJE-1392 15.0 cm, serrated



FORCEPS / THUMB

Gillies Dissecting



DJE-1393 15.0 cm, cross serrated, 1x2 teeth

General



DJE-1394 15.0 cm, cross serrated, TC

Cushing



DJE-1395 17.0 cm, cross serrated, TC

Cushing



DJE-1396 17.0 cm, serrated, Titanium



FORCEPS / THUMB

Cushing



DJE-1397 17.0 cm, angled,
cross serrated, TC

Taylor Cushing



DJE-1398 18.5 cm,
cross serrated,
bayonet, TC

Bonney



DJE-1399 17.5 cm, cross serrated,
1x2 teeth, TC

Potts-Smith



DJE-1400 18.0 cm, cross serrated, TC
DJE-1401 21.0 cm, cross serrated, TC
DJE-1402 25.0 cm, cross serrated, TC



FORCEPS / THUMB

Gerald



- DJE-1403 18.0 cm, cross serrated, 1x2 teeth, TC
- DJE-1404 18.0 cm, cross serrated, TC

Skinner



- DJE-1405 23.0 cm, cross serrated, 1x2 teeth, TC

DeBakey Titanium



- DJE-1406 15.0 cm, 1.5 mm tip, Titanium
- DJE-1407 19.5 cm, 1.5 mm tip, Titanium
- DJE-1408 24.0 cm, 1.5 mm tip, Titanium
- DJE-1409 15.0 cm, 2.3 mm tip, Titanium
- DJE-1410 19.5 cm, 2.3 mm tip, Titanium
- DJE-1411 24.0 cm, 2.3 mm tip, Titanium

DeBakey



- DJE-1412 16.0 cm, 1.5 mm jaw
- DJE-1413 20.0 cm, 1.5 mm jaw
- DJE-1414 24.0 cm, 1.5 mm jaw
- DJE-1415 16.0 cm, 2.0 mm jaw
- DJE-1416 20.0 cm, 2.0 mm jaw
- DJE-1417 24.0 cm, 2.0 mm jaw
- DJE-1418 30.0 cm, 2.0 mm jaw
- DJE-1419 16.0 cm, 2.7 mm jaw
- DJE-1420 20.0 cm, 2.7 mm jaw
- DJE-1421 24.0 cm, 2.7 mm jaw
- DJE-1422 30.0 cm, 2.7 mm jaw
- DJE-1423 16.0 cm, 3.5 mm jaw
- DJE-1424 20.0 cm, 3.5 mm jaw
- DJE-1425 24.0 cm, 3.5 mm jaw
- DJE-1426 30.0 cm, 3.5 mm jaw



FORCEPS / THUMB

Lucae-Jansen



DJE-1427 16.0 cm, serrated

Jansen Gruenwald



DJE-1428 20.0 cm, serrated

Lucae



DJE-1429 14.0 cm, serrated











FORCEPS / THUMB

Jewelers diamond dusted



- DJE-1430  11.5 cm, # 5A, very fine pointed tips
- DJE-1431  11.5 cm, # 3, straight, delicate tips
- DJE-1432  11.5 cm, # 4, straight, very fine tips
- DJE-1433  11.0 cm, # 5, straight, ultra fine tips
- DJE-1435  11.0 cm, # 5/45, angled, ultra fine tips
- DJE-1436  11.5 cm, # 7, curved, delicate tips

Jewelers

- DJE-1437  11.5 cm, # 1, straight, standard tips
- DJE-1438  11.5 cm, # 3, straight, delicate tips
- DJE-1439  11.0 cm, # 5A, very fine pointed tips
- DJE-1440  10.5 cm, # 3C, straight, delicate tips
- DJE-1441  11.0 cm, # 4, straight, very fine tips
- DJE-1442  11.0 cm, # 5, straight, ultra fine tips
- DJE-1443  11.0 cm, # 5/45, angled, ultra fine tips
- DJE-1445  11.0 cm, # 7, curved, delicate tips



Micro-Forceps

. round handle, with tying platform



- DJE-1446 15.0 cm, straight, smooth
- DJE-1447 18.0 cm, straight, smooth

- DJE-1448 15.0 cm, curved, smooth
- DJE-1449 18.0 cm, curved, smooth



MICRO SURGERY INSTRUMENTS

Vessel Dilator Angled

- . 11cm/4¼" length
- . Tips 0.2mm Ø



DJE-1450

Vessel Dilator Angled

- . 11cm/4¼" length
- . Tips 0.3mm Ø



DJE-1451

Vessel Dilator

- . 13.5cm/5¾" length
- . Tips 0.2mm Ø



DJE-1452

Vessel Dilator

- . 13.5cm/5¾" length
- . Curved tips 0.2mm Ø



DJE-1453



MICRO SURGERY INSTRUMENTS

Vessel Dilator Long


- . 18cm/7" length
- . Tips 0.3mm \emptyset



DJE-1454

Vessel Dilator Long

- . 18cm/7" length
- . Curved tips 0.2mm \emptyset



DJE-1455

Round Bodied Vessel Dilator

- . 15cm/6" length
- . Straight tips 0.3mm \emptyset
- . Balanced forceps



DJE-1456

Round Bodied Vessel Dilator

- . 15cm/6" length
- . Curved tips 0.2mm \emptyset
- . Balanced forceps



DJE-1457



MICRO SURGERY INSTRUMENTS

Vessel Dilator Long

- . 18cm/7" length
- . Straight tips 0.3mm Ø
- . Balanced forceps



DJE-1458

Round Bodied Vessel Dilator

- . 18cm/7" length
- . Curved tips 0.3mm Ø
- . Balanced forceps



DJE-1459

Pierse Forceps

- . 12cm/4¾" length
- . Tips 0.3mm



DJE-1460

Micro Adson Fixation Forceps

- . 12cm/4¾" length
- . 1 x 2 micro teeth
- . Tips 0.3mm



DJE-1461



MICRO SURGERY INSTRUMENTS

Round Bodied Forceps

- . 15cm/6" length
- . Straight, with platform
- . Tips 0.3mm \emptyset
- . Balanced forceps



DJE-1462

Round Bodied Forceps

- . 15cm/6" length
- . Curved, with platform
- . Tips 0.3mm \emptyset
- . Balanced forceps



DJE-1463

Round Bodied Angled Forceps

- . 15cm/6" length
- . Angled tips 0.2mm \emptyset
- . Balanced forceps



DJE-1464



MICRO SURGERY INSTRUMENTS

Round Bodied Forceps

- . 18cm/7" length
- . Straight, with platform
- . Tips 0.5mm Ø
- . Balanced forceps



DJE-1465

Round Bodied Forceps

- . 18cm/7" length
- . Curved, with platform
- . Tips 0.5mm Ø
- . Balanced forceps



DJE-1466

Round Bodied Angled Forceps

- . 18cm/7" length
- . Angled tips 0.2mm Ø
- . Balanced forceps



DJE-1467




MICRO SURGERY INSTRUMENTS

Straight Scissors Flat Handle

- . 15cm/6" length
- . Blades 19mm long



Adventitia Scissors with pointed tips



DJE-1469

Available with Micro Serrations for nerve cutting - **Straight**


DJE-1468

Curved Scissors Flat Handle

- . 15cm/6" length
- . Blades 19mm long



Dissecting Scissors with rounded tips



DJE-1471

Available with Micro Serrations for nerve cutting - **Curved**


DJE-1470

Straight Scissors Flat Handle

- . 15cm/6" length
- . Blades 15mm long



Adventitia Scissors with pointed tips




DJE-1472

Curved Scissors Flat Handle

- . 15cm/6" length
- . Blades 15mm long



Dissecting Scissors with rounded tips



DJE-1473



MICRO SURGERY INSTRUMENTS

Straight Scissors Round Handle

- . 15cm/6" length
- . Blades 12mm long


DJE-1474

- . 18cm/7" length
- . Blades 14mm long



DJE-1475

Adventitia scissors with pointed tips



Curved Scissors Round Handle

- . 15cm/6" length
- . Blades 12mm long


DJE-1476

- . 18cm/7" length
- . Blades 14mm long



DJE-1477

Dissecting scissors with rounded tips



Blondeel Scissors Round Handle

- . Blades 15mm long
- . Curved blades



DJE-1478 15cm/6" length
 DJE-1479 18cm/7" length

These Scissors were designed by professor Phillip Blondeel for use in breast reconstruction utilising a perforator flap. A combination of ring handle design for preparation of fine structures, the round handle micro instrument handle allows the user to perform fine micro surgery.



MICRO SURGERY INSTRUMENTS

Barraquer Micro Needle Holder

- . Round bodied
- . 0.4mm jaws



Without Lock	Straight	Curved
13cm/5"	DJE-1480	DJE-1481
15cm/6"	DJE-1482	DJE-1483
18cm/7"	DJE-1484	DJE-1485

Castroviejo Micro Needle Holder

- . Flat handled
- . 0.4mm jaws



Without Lock	Curved
14cm/5½"	DJE-1486
18cm/7"	DJE-1487

Double Action Micro Needle Holder

- . Round bodied
- . 0.4mm jaws




Without Lock	Curved
15cm/6"	DJE-1488
18cm/7"	DJE-1489



ULTRA FINE MICRO INSTRUMENTS

Ultra Fine Jewellers Forceps

- . 11cm/4¼" length
- . Tips 0.1mm Ø



DJE-1490

Ultra Fine Jewellers Forceps


- . 13.5cm/5¾" length
- . Tips 0.1mm Ø



DJE-1491

Ultra Fine Round Bodied Forceps


- . 15cm/6" length
- . Straight tips 0.1mm Ø
- .. Straight tips 0.1mm Ø



DJE-1492

Ultra Fine Round Bodied Forceps

- . 18cm/7" length
- . Straight tips 0.1mm Ø
- . Balanced forceps



DJE-1493



ULTRA FINE MICRO INSTRUMENTS

Ultra Fine Vessel Dilator Angled


- . 11cm/4¼" length
- . Tips 0.1mm Ø



DJE-1494

Ultra Fine Vessel Dilator Angled


- . 13.5cm/5¾" length
- . Curved tips 0.2mm Ø



DJE-1495

Ultra Fine Round Bodied Vessel Dilator


- . 15cm/6" length
- . Straight tips 0.1mm Ø
- . Balanced forceps



DJE-1496

Ultra Fine Round Bodied Vessel Dilator

- . 15cm/6" length
- . Curved tips 0.2mm Ø
- . Balanced forceps



DJE-1497



ULTRA FINE MICRO INSTRUMENTS

Ultra Fine Straight Scissors Flat Handle

- . 15cm/6" length
- . Blades 14mm long



Adventitia scissors with pointed tips



DJE-1498

Ultra Fine Curved Scissors Flat Handle

- . 15cm/6" length
- . Blades 12mm long



Dissecting scissors with rounded tips



DJE-1499

Ultra Fine Straight Scissors Round Handle

- . Blades 14mm long



Dissecting scissors with pointed tips

DJE-1500 15cm/6" length

DJE-1501 18cm/7" length

Ultra Fine Curved Scissors Round Handle

- . Blades 12mm long



Dissecting scissors with rounded tips

DJE-1502 15cm/6" length

DJE-1503 18cm/7" length



ULTRA FINE MICRO INSTRUMENTS

Ultra Fine Curved Micro Needle Holder

- . 14cm/5½" length
- . Flat handled
- . 0.2mm jaws



DJE-1504

Ultra Fine Curved Micro Needle Holder

- . Round bodied
- . 0.2mm jaws



DJE-1505 15cm/6" length

DJE-1506 18cm/7" length

Ultra Fine Set Complete



DJE-1507



DISSECTORS

Freers Elevator/Dissectort

. 18.5cm/7¼" length



DJE-1508

Watson Cheyne Dissector



	Size
DJE-1509	18cm/7"
DJE-1510	23cm/9"
DJE-1511	25.5cm/10"

MacDonald Elevator/Dissector



	Size
DJE-1512	19cm/7½"
DJE-1513	23cm/9"

Howarth Elevator

. 21.5cm/8½" length



DJE-1514

Mitchels Trimmer

. 16cm/6¼" length



DJE-1515



RETRACTORS

Ferreira retractor

DJE-1516 11.0 cm x 2.2 cm, face lift retractor with fiberoptic



Flap retractor

DJE-1517 11.0 cm, 2.2 cm blade



Tebbetts Style

- . T-Mammoplasty retractor with smoke
- . evacuation tube and fiberoptical illumination,
- . with Storz, HSW, ACMI adapter included

- DJE-1518 16 x 80 mm blade, without teeth
- DJE-1519 24 x 90 mm blade, without teeth
- DJE-1520 30 x 90 mm blade, without teeth
- DJE-1521 24 x 100 mm blade, without teeth
- DJE-1522 25 x 120 mm blade, without teeth
- DJE-1523 30 x 120 mm blade, without teeth
- DJE-1524 25 x 150 mm blade, without teeth
- DJE-1525 30 x 150 mm blade, without teeth
- DJE-1526 36 x 150 mm blade, without teeth
- DJE-1527 40 x 170 mm blade, without teeth
- DJE-1528 30 x 190 mm blade, without teeth

- DJE-1529 16 x 80 mm blade, with teeth
- DJE-1530 24 x 90 mm blade, with teeth
- DJE-1531 30 x 90 mm blade, with teeth
- DJE-1532 24 x 100 mm blade, with teeth
- DJE-1533 25 x 120 mm blade, with teeth
- DJE-1534 30 x 120 mm blade, with teeth
- DJE-1535 25 x 150 mm blade, with teeth
- DJE-1536 30 x 150 mm blade, with teeth
- DJE-1537 36 x 150 mm blade, with teeth
- DJE-1538 40 x 170 mm blade, with teeth
- DJE-1539 30 x 190 mm blade, with teeth





RETRACTORS

Deaver



DJE-1540 20.0 cm, retractor with fiberoptic

Freeman Flap



DJE-1541 20.5 cm,
retractor, 12 mm blade width
DJE-1542 20.5 cm,
retractor, 25 mm blade width

Multiple approach retractor



DJE-1543 19.0 cm, 23 mm blade



RETRACTORS

Biggs Mammoplasty Retractor

- DJE-1544 standard blade, 50 mm
- DJE-1545 standard blade with fiber-optics, 50 mm
- DJE-1546 standard blade with fiber-optics and suction tube, 50 mm
- DJE-1547 narrow blade, 31 mm
- DJE-1548 narrow blade with fiber-optics, 31 mm
- DJE-1549 narrow blade with fiber-optics and suction tube, 31 mm



Converse

- DJE-1550 11.0 cm, nasal retractor, 11x44 mm blade
- DJE-1551 11.0 cm, nasal retractor, 11x44 mm blade, with fiber optics



Aufricht

- DJE-1552 16.0 cm, nasal retractor, 45 mm blade
- DJE-1553 16.0 cm, nasal retractor, 45 mm blade, with fiber optics



Cottle

- DJE-1554 20.0 cm, retractor with knife guard





RETRACTORS

Weitlaner



- DJE-1555 11.5 cm, retractor, 2x3 sharp prongs
- DJE-1556 14.0 cm, retractor, 3x4 sharp prongs

Swivel double hook



- DJE-1557 8.0 cm, sharp, 10 mm

Sheen



- DJE-1558 16.0 cm, nasal retractor

Spatula

. semi-malleable



- DJE-1559 31.5 cm, 2.5x15.5 cm blade



RETRACTORS

Hamra



DJE-1560 16.0 cm, double ended retractor

Cottle



DJE-1561 14.0 cm, retractor, double hook, sharp

Fomon



DJE-1562 16.0 cm, alar retractor, double hook with ball ends

Converse Alar



DJE-1563 10.0 cm, retractor double ended, 10mm / 14mm, blunt

Martinez Kilner



DJE-1564 11.5 cm, retractor double hook with thimble, 10mm wide, sharp



RETRACTORS

U.S. Army



DJE-1565 21.0 cm, retractor,
double ended, pair

Maxwell Flap



DJE-1566 10.0 cm x 22 mm, four prongs retractor
DJE-1567 10.0 cm x 38 mm, four prongs retractor

Freeman Rake



DJE-1568 18.0 cm x 22 mm,
four prongs retractor
DJE-1569 18.0 cm x 38 mm,
four prongs retractor

Anderson Bear Claw



DJE-1570 11.0 cm, 3 prongs, sharp
DJE-1571 11.0 cm, 5 prongs, sharp



RETRACTORS

Freeman



DJE-1572 10.0 cm, double skin hook, sharp

Yancoskie lift retractor



DJE-1573 15.0 cm, 6 prongs

Lift retractor



DJE-1574 18.0 cm, 7 prongs



RETRACTORS

Gorney Rake



DJE-1575 7.5 cm, retractor with spring, 38mm wide

Rollet



DJE-1576 13.0 cm, retractor, delicate,
4 sharp prongs

Miller Rake



DJE-1577 13.5 cm, tissue retractor, delicate,
9 sharp prongs

Ragnell



DJE-1578 15.0 cm, double ended retractor,
delicate

Ragnell-Davis



DJE-1579 14.5 cm,
double ended retractor



RETRACTORS

Senn



- DJE-1580 17.0 cm, double ended retractor, sharp
- DJE-1581 17.0 cm, double ended retractor, blunt


Holtzheimer

. self retaining retractor, 4 teeth



- DJE-1582 10.0 cm, blunt, straight
- DJE-1583 10.0 cm, sharp, straight
- DJE-1584 10.0 cm, blunt, curved
- DJE-1585 10.0 cm, sharp, curved

Farr Spring Retractor



- DJE-1586 14.0 cm x 2.0 cm
- DJE-1587 14.0 cm x 3.0 cm



RETRACTORS

Kasden



DJE-1588 16.5 cm,
double ended retractor

Kawamoto



DJE-1589 16.0 cm,
double ended retractor,
insulated, 8mm/10mm blades

Kawamoto




DJE-1590 16.0 cm,
double ended retractor,
insulated, 14mm/16mm blades

Wise



DJE-1591 16.5 cm, orbit retractor,
double ended

Hajek's



DJE-1592 18.0 cm, cheek retractor



RETRACTORS

Sheen



DJE-1593 15.0 cm, retractor

Caroline



DJE-1594 11.0 cm, finger retractor

Cronin



DJE-1595 15.5 cm, lip and nostril retractor

Cronin



DJE-1596 16.5 cm, cheek retractor

Fischl



DJE-1597 16.0 cm,
"Shoe Horn" retractor



SELF RETAINING RETRACTORS

Alms Retractor

- . 7cm/2¾" length
- . 4 x 4 teeth blunt
- . 4mm depth prongs

DJE-1598



Baby Weitlaner S/R Retractor

- . 11.5cm/4½" length
- . 2 x 3 blunt prongs

DJE-1599



Wests S/R Retractor

- . 14cm/5½" length
- . 3 x 4 sharp prongs

DJE-1600 *



Weitlaner S/R Retractor

- . 14cm/5½" length
- . 3 x 4 blunt prongs

DJE-1601

Mollison S/R Retractor - Child

- . 13cm/5¼" length
- . 2 x 2 sharp prongs 12mm depth
- . Curved

DJE-1602



Mollison S/R Retractor - Adult

- . 17cm/6¾" length
- . 4 x 4 sharp prongs 14mm depth
- . Curved

DJE-1603





RETRACTORS

Kilner-Miniature Pattern Catspaw Retractor

. 15cm/6" length



DJE-1604

Senn-Miller Double Ended Retractor

. 16cm/6¼" length



DJE-1605 Sharp prongs
DJE-1606 Blunt prongs

Ragnell-Davies Retractor

. 14cm/5½" length



DJE-1607

Mini-Langenbeck Retractor

. 16cm/6¼" length



DJE-1608 10 x 6mm
DJE-1609 17 x 5mm

Desmarres Saddle Retractor

. 16cm/6¼" length
. 12mm width



DJE-1610

Langenbeck Retractor

. 21.5cm/8½" length

Code	Blade Depth	Size Width
DJE-1611	22mm x 6mm	Small
DJE-1612	45mm x 13mm	Medium
DJE-1613	65mm x 25mm	Large





RETRACTORS & SKIN HOOKS

Kilner Detachable Skin Hook

. 15cm/6" length



DJE-1614

Gillies Skin Hook

. 18cm/7" length



DJE-1615 Small size

DJE-1616 Large size

Klienert-Kutz Skin Hook

. 16cm/6¼" length



DJE-1617 Small size

DJE-1618 Large size

Joseph Double Skin Hook

. 16cm/6¼"



DJE-1619

Rollet/Blair Mini-Catspaw Retractor

. 13.5cm/5" length
. Semi-sharp prongs



DJE-1620

Catspaw Retractors

. 16cm/6¼" length
. Semi-sharp prongs



DJE-1621



MAMMOPLAST BREAST RETRACTOR

- . Self retaining ratchet
- . Adjustable to hold breast during surgery




DJE-1622



BREAST RETRACTORS

Tebbetts Fibre Optic Breast Retractor

. With fibre optic illumination and suction



- . Fibre optic illumination with suction for smoke evacuation
- . Angled Retractor tip with blunt teeth
- . Polished mirror finish on blade to reflect light inside the pocket



SCALPEL HANDLES & MARKING PENS

Barrons Scalpel Blade Handle No.3

. 13cm/5¼" length

DJE-1631



Beaver Blade Holder

. 13cm/5¼" length

DJE-1632



Eckoff Mapping Pen

. 14.5cm/5¾" length

DJE-1633



Sommerlad Marking Pen

. 15cm/6" length

DJE-1634



Castroviejo Measuring Caliper

. 7cm/2¾" length
. Measuring range 0 to 25mm



DJE-1635

Ink Pot

. Stainless steel



DJE-1636

Castroviejo Measuring Caliper

. 7cm/2¾" length
. Measuring range 0 to 25mm



DJE-1637



SKIN GRAFT KNIVES

Watson Skin Graft Knife

. 30.5cm/12" length



Right handed

DJE-1638


Left handed

DJE-1639

The roller bar is smooth and does not oscillate or rotate in use. The graft thickness is controlled by slackening the locking nut and rotating the adjuster nut until the required thickness and then set by tightening the locking nut. The graft thickness is consistent as the roller bar can only be adjusted in parallel.

Braithwaite Skin Graft Knife

. 30.5cm/12" length



Right handed

DJE-1640

Left handed

DJE-1641

The grooved roller bar oscillates from side to side as well as rotating. The graft thickness is controlled by rotating the adjuster nuts at the end of knife, this will move the roller bar up or down.

Cobbett Skin Graft Knife

. 30.5cm/12" length



Right handed

DJE-1642


Left handed

DJE-1643

The grooved roller bar only oscillates from side to side. The adjustment for the thickness of the graft is the same as the Braithwaite knife.

Humby Skin Graft Knife

. 36.5cm/14 1/4" length



Right handed

DJE-1644

Left handed

DJE-1645

The grooved roller bar rotates and oscillates from side to side. The graft thickness is controlled by adjusting the thumb screws on top of the support pillars.

Silver Skin Graft Knife

. 19cm/7 1/2" length



Right handed

DJE-1646

Knife

Cobbett

Watson

Braithwaite

Silver

Left handed

DJE-1647

Blades code

DJE-1648

DJE-1649

DJE-1650

DJE-1651

Tufnol Skin Graft Board

200mm x 100mm/ 8 x 4"
Bevelled edge



DJE-1652



AREOLA MARKERS

Freeman Areola Markert

Code	Diameter
DJE-1653	36mm Ø
DJE-1654	38mm Ø
DJE-1655	40mm Ø
DJE-1656	42mm Ø
DJE-1657	45mm Ø



Freeman Double Ended Areola Marker

Code	Diameter
DJE-1658	38/42mm Ø
DJE-1659	42/45mm Ø



Grossman Areola Marker

- . 4 sizes in 1
- . Marking 40mm/44mm/48mm & 52mm Ø



DJE-1660

Galvao Nipple Marker

Code	Diameter
DJE-1661	40mm Ø standard
DJE-1662	50mm Ø large



McKissock Keyhole Breast Reduction Marker

Code	Diameter
DJE-1663	38mm Ø
DJE-1664	42mm Ø
DJE-1665	45mm Ø



- . Precisely marks amount of Tissue for removal
- . Used together with Freeman Marker
- . Arms are graduated in centimeters for marking tissue to be removed



FIBRE OPTIC LIGHT LEAD

Fibre Optic Light Lead



Code	Length	Diameter
DJE-1666	1800mm	3.5mm Ø
DJE-1667	2300mm	3.5mm Ø
DJE-1668	3000mm	3.5mm Ø
DJE-1669	1800mm	4.8mm Ø
DJE-1670	2300mm	4.8mm Ø
DJE-1671	3000mm	4.8mm Ø



ADAPTORS FOR LIGHT LEAD

Lightsource Fittings



Wolf DJE-1672



ACMI British Standard DJE-1673



Olympus - Standard DJE-1674



Mercian/Storz DJE-1675



Olympus DJE-1676

Retractor Fitting



DJE-1677