

A Good Helper of Vision Screening

Chrops vision screener can screening shortsightedness,hyperopia, astigmatism,PD,pupil size and anisometropia fast and accurately. And it provides reference for prevention activities.


Test Both Eyes at One Time

Refractive diopter, PD, PD diameter, fixation asymmetry, anisometropia and other data of both eyes can be received quickly with just one test.




Rapid Measurement

1m test distance,non-contact test, safe and more sanitary,smart far and near focus prompts.The fastest measurement just need 0.1sec.



Data Transmission

Data can be bulk imported and exported by connecting with USB stick, USB cable,WIFI.




Bluetooth Print

The external bluetooth printer is small but sensitive. If there is no need to print out, also can get the measurement result by scanning Wechat QR code.




Built-in Rechargeable Battery

Built-in rechargeable battery with high power capacity; adjustable brightness; 3 levels of settings for intelligent screen saver time; auto power-off if not using for a long time.




Wide Range of Measurement

The range of measurement is from 6 months to 100-year-old. In the back, there provides dynamic led and music to attract kids' eyes for easy measurements.



Digital List, Easy & Quick

You can do paging query also quick search by names. Testing reports can be created after connecting with matching PC softwares.



Small Size, Light Weight, Easy to Carry