

Specification

Acquisition Modes	Non-mydrriatic/mydrriatic Anterior photography (FFA) / (FAF)
Field View	50°
Working Distance	35mm
Minimum Pupil	≥3.3mm
Focus Modes	Manual/Auto
Alignment Modes	Double dots auxiliary
Exposure	Manual /Auto
Photography	SLR Camera
Image definition	24 Megapixel
Compensation	±25D
Base Movement	Front/Back : 60mm ; Right/Left : 100mm ; Up/Down : 30mm
Left/Right angle	Left/Right 30°
Pitch Angle	up : 15° , down : 12°
Fixation	External /Internal (any position available)
Screen	10.1 inch touch screen
DICOM 3.0	Support
Powerful software	Image process ;Auto mosaic ;Follow-up ;LAN, FTP transport ;
Size	370mm (L) ×520mm (W) ×600mm (H)
Weight	35kg
Power	100 ~ 240V 50/60Hz 200VA

RetiCam 3000 Plus Multifunction Fundus Camera



Technical Feature :

- * Non-mydrriatic, automatic focus and exposure.
- * 10.1 inch touch screen, easy operation.
- * Compact design with smaller computer host.
- * Red-free, color, FFA and FAF are available.
- * Powerful image process system.
- * Expansion port include DICOM 、 LAN 、 USB, build remote and AI platform easily.
- * 2400Mega Pixel , present minor lesions .

