



Portable Slit Lamp

Introduction



Technical Parameter

Type Binocular Portable Slit Lamp	Angle of Convergence 13°
Magnification Model LS-1A 10X	Model LS-1B 10X 16X
Working Distance 80mm (16X) 100mm(10X)	Pupil Adjustment Range 49-75 mm
Diopter Adjustment ±7D	Slit Angle ±60°
Filter Neutral Density, Red-free, Cobalt Blue	Light Resource Warm LED, Color Temperature 3500k
Illumination Slit Width 0.1mm, 0.2mm, 0.8mm & Φ1mm, Φ5mm, Φ12mm circle	
Slit Length 12mm(with width 0.1mm, 0.8mm), 8mm(with width 0.2mm)	
Illumination Control Continuously Variable from Low to Full Brightness(≥250 000 lux)	
Power (Battery)	
LED Working Voltage AA, 15VA	Battery Type AA Dry Cell / Rechargeable Batteries
Continuous Operating Time ≥ 4 hours	
Dimensions and Weight Dimensions 320*310*205 mm Net Weight Approx. 835g	
Standard Accessories: Carrying Case, Battery Charger, Forehead Support, Instrument Manual	