

CS8.0



Patent No.: 2019305650546

European Patent No.: 007081682-0001

USA Patent No.: US D900,416 S



SPECIFICATION

Tank Capacity: 8000ml

Maximum: 7000ml

Power Supply: 450W AC100-120V 50/60Hz

450W AC220-240V 50/60Hz

Digital Timer Settings: 5min-10min-15min-20min-25min-30min

Frequency: 46KHz \pm 3KHz

Tank Material: Stainless steel SUS304

Tank Size: 33.1x23.8x11.6 cm (13.0"x9.4"x4.6")

Unit Size: 46.6x35x22.9 cm (18.3"x13.8"x9.0")

Inner Carton Size: 52.5x39x27.5 cm (20.7"x15.3"x10.8")

Master Carton Size: 54x40.5x30.5 cm (21.3"x16.0"x12.0")

Net Weight: 14.3kg(31.5lb)

Unit Weight of Inner Carton: 16.0kg(35.2lb)

Qty per Master Carton: 1pc/ctn

Gross Weight: 17.1kg(37.7lb)